

C-PAK™ CMC DISSOLVABLE NASAL DRESSING



Affordable and easy to use

C-Pak™ CMC Dissolvable Nasal Dressing is composed of crosslinked CMC, which is hydrated to a gel immediately prior to use. The gel is eliminated within 14 days with adequate hydration, or can be suctioned or irrigated from the cavity earlier at the physicians discretion.

CMC (carboxymethylcellulose) is a plant derivative, with high viscosity, non-toxic, non-allergenic properties.



Reconstituted CMC hydrogel

ADVANTAGES OF C-PAK™ DRESSING

- 14 cc of hydrogel
- Innovative mixing technique
- Absorbs 16 times its weight
- Hydrates with sterile saline
- Flexible cannula

ordering information

qty

1524000

C-Pak™ CMC Dissolvable Nasal Dressing

1

Medtronic ENT

Medtronic USA, Inc.

6743 Southpoint Drive North
Jacksonville, FL 32216
USA

www.MedtronicENT.com

Telephone: (904) 296-9600

Fax: (904) 296-9666

International Telephone Numbers

Australia 1-800-668-670

Belgium 32-2456-09-00

Canada 1-800-217-1617

China 86-21-50800998

France 33-470-679-800

Germany 49-211-5293-209

Hong Kong 852-2919-1312

India 91-22-26836733

Italy 39-02-24137-324

Japan 81-6-4795-1506

Korea: 82-2-3404-3600

Lebanon: 961-1-370-670

Luxembourg: 32-2456-09-09

Netherlands: 31-45-566-8000

Poland: 48-22-465-6942

Singapore: 65-6776-6255

Spain: 34-91-625-05-40

UK: 44-1923-2-212213

USA: 1-904-296-9600

For further information,
please call Medtronic ENT at
(904) 296-9600. You may also
consult our website
at www.MedtronicENT.com.

This literature is intended for the exclusive use of physicians. Rx only.

©2008, Medtronic, Inc. C-Pak™ is a registered trademark of Medtronic, Inc.

Sinu-Foam™ is a registered trademark of ArthroCare Corporation.

890916EN 05.08 2008-413



Medtronic

Alleviating Pain · Restoring Health · Extending Life